



Ablaze Glass Works Pvt. Ltd.



Quick Setup Lab Reactor


*Interchangeability & Scalability
Reaction Station for study & Analysis*

250 to 5000mL



www.ablazeglassworks.com
marketing@ablazeglassworks.com

Your Quick Setup Reaction Station



Jacketed Condenser

5 Neck DN100 Lid
Common for all reactor

Double Wall/Triple
Reactor for Desired Reaction
250 mL~5000 mL

Quick Release Coupling
Reactor Interchangeability

Hose Pipe Connectivity
Avoid Stress on reactor

Manifold for Connectivity
with Utility

**Single Station
for Many Reactions.
Scale your Reactions**

Pressure Equalizing
Additional Funnel

Choice of OHS


Magnetic Coupling
upto 1800 RPM

SS 304 Scaffold

PTFE sleeved SS316 Stirrer
with Choice of Blade

Zero Hold Drain Valve

Automation



Quick Setup Reaction Station

Choice of Reaction Vessel
250 mL to 5000 mL



Jacketed
Condenser



Quick Release Coupling
for Easy Swapping



Additional
Funnel

Choice of
Stirrer Blade



Zero Hold Flush
Bottom Valve





Rodavis Joint



Magnetic Coupling for Higher RPM



PTFE Baffle Available



Bola Connectors



Quick Setup Lab Reactor




Overheated Stirrers



Compatible with All International Brands



Automation :
On Request





Lab Products

High Pressure Glass Reactor <<

Lab Rotary Evaporator <<

Lab Gas Scrubber Unit <<

Lab Short Path Distillation <<

Glass Schlenk Manifold <<



Industrial Products

Industrial Components <<

Industrial Units <<

Made to Design Components <<

Contract Manufacturing <<

Ablaze Glass Works Pvt. Ltd.



E-52, Sardar Industrial Estate,
Vadodara, Gujarat, India.



www.ablazeglassworks.com



marketing@ablazeglassworks.com